


Product datasheet for **GA105715**

Human SKAP55 (SKAP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SKAP1 CRISPRa kit - CRISPR gene activation of human src kinase associated phosphoprotein 1
Symbol:	SKAP1
Components:	GA105715G1, SKAP55 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105715G2, SKAP55 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105715G3, SKAP55 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86WV1
Synonyms:	HEL-S-81p; SCAP1; SKAP55
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)