


Product datasheet for **GA105716**

Human DNAH17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNAH17 CRISPRa kit - CRISPR gene activation of human dynein axonemal heavy chain 17
Symbol:	DNAH17
Components:	GA105716G1, DNAH17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105716G2, DNAH17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105716G3, DNAH17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UFH2</u>
Synonyms:	axonemal dynein heavy chain; DNAHL1; DNEL2; dynein, axonemal, heavy chain 17; dynein, axonemal, heavy chain like 1; dynein, axonemal, heavy like 1; dynein, axonemal, heavy polypeptide 17; FLJ40457; MGC132767; MGC138489
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)