


## Product datasheet for **GA105728**

### Human KCNK5 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNK5 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 5
Symbol:	KCNK5
UniProt ID:	<u><a href="#">O95279</a></u>
Synonyms:	K2p5.1; KCNK5b; TASK-2; TASK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)