

Product datasheet for **GA105737**

Human **DYNLL1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DYNLL1 CRISPRa kit - CRISPR gene activation of human dynein light chain LC8-type 1
Symbol:	DYNLL1
UniProt ID:	P63167
Synonyms:	DLC1; DLC8; DNCL1; DNCLC1; hdlc1; LC8; LC8a; PIN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)