

Product datasheet for **GA105741**

Human EIF3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3A CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit A
Symbol:	EIF3A
UniProt ID:	Q14152
Synonyms:	EIF3; eIF3-p170; eIF3-theta; EIF3S10; P167; p180; p185; TIF32
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)