

Product datasheet for **GA105743**

Human eIF3s8 (EIF3C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3C CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit C
Symbol:	EIF3C
UniProt ID:	Q99613
Synonyms:	eIF3-p110; EIF3CL; EIF3S8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)