

## Product datasheet for **GA105744**

### Human EIF3D activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | EIF3D CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit D |
| Symbol:       | EIF3D  |
| UniProt ID:   | <u><a href="#">O15371</a></u>  |
| Synonyms:     | eIF3-p66; eIF3-zeta; EIF3S7  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                  |

#### Product images:



[View online >](#)