

Product datasheet for **GA105745**

Human EIF3F activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3F CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit F
Symbol:	EIF3F
UniProt ID:	O00303
Synonyms:	eIF3-p47; EIF3S5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)