

Product datasheet for **GA105746**

Human EIF3S4 (EIF3G) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3G CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit G
Symbol:	EIF3G
UniProt ID:	<u>O75821</u>
Synonyms:	eIF3-delta; EIF3-P42; eIF3-p44; EIF3S4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)