


Product datasheet for **GA105747**

Human EIF3S3 (EIF3H) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3H CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit H
Symbol:	EIF3H
UniProt ID:	<u>O15372</u>
Synonyms:	eIF3-gamma; eIF3-p40; EIF3S3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)