

Product datasheet for **GA105748**

Human EIF3S2 (EIF3I) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3I CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit I
Symbol:	EIF3I
Components:	GA105748G1, EIF3S2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105748G2, EIF3S2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105748G3, EIF3S2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q13347</u>
Synonyms:	eIF3-beta; eIF3-p36; EIF3S2; PRO2242; TRIP-1; TRIP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)