

Product datasheet for **GA105749**

Human EIF3S1 (EIF3J) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3J CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit J
Symbol:	EIF3J
Components:	GA105749G1, EIF3S1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105749G2, EIF3S1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105749G3, EIF3S1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q75822</u>
Synonyms:	eIF3-alpha; eIF3-p35; EIF3S1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

