

## Product datasheet for **GA105751**

### Human eIF4GII (EIF4G3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	EIF4G3 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 4 gamma 3
Symbol:	EIF4G3
UniProt ID:	<a href="#">O43432</a>
Synonyms:	eIF-4G 3; eIF4G 3; eIF4GII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)