

## Product datasheet for **GA105785**

### Human DR3 (TNFRSF25) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF25 CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 25
Symbol:	TNFRSF25
UniProt ID:	<a href="#">Q93038</a>
Synonyms:	APO-3; DDR3; DR3; LARD; TNFRSF12; TR3; TRAMP; WSL-1; WSL-LR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)