

Product datasheet for **GA105806**

Human TWEAK (TNFSF12) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFSF12 CRISPRa kit - CRISPR gene activation of human TNF superfamily member 12
Symbol:	TNFSF12
UniProt ID:	O43508
Synonyms:	APO3L; DR3LG; TWEAK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)