

Product datasheet for **GA105822**

Human Dcr3 (TNFRSF6B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF6B CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 6b
Symbol:	TNFRSF6B
UniProt ID:	<u>O95407</u>
Synonyms:	DCR3; DJ583P15.1.1; M68; M68E; TR6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)