

Product datasheet for **GA105831**

Human TNFRSF18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF18 CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 18
Symbol:	TNFRSF18
UniProt ID:	Q9Y5U5
Synonyms:	AITR; CD357; GITR; GITR-D
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)