

Product datasheet for **GA105840**

Human DcR1 (TNFRSF10C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF10C CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 10c
Symbol:	TNFRSF10C
UniProt ID:	O14798
Synonyms:	CD263; DCR1; DCR1-TNFR; LIT; TRAIL-R3; TRAILR3; TRID
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)