

Product datasheet for **GA105863**

Human SAP30 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SAP30 CRISPRa kit - CRISPR gene activation of human Sin3A associated protein 30
Symbol:	SAP30
UniProt ID:	O75446
Synonyms:	30kD; histone deacetylase complex subunit SAP30; Sin3 corepressor complex subunit SAP30; Sin3-associated polypeptide, 30kDa; Sin3A-associated protein, 30kDa
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)