

## Product datasheet for **GA105876**

### Human TMEM11 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | TMEM11 CRISPRa kit - CRISPR gene activation of human transmembrane protein 11 |
| Symbol:       | TMEM11  |
| UniProt ID:   | <a href="#">P17152</a>  |
| Synonyms:     | C17orf35; PM1; PMI  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)       |

#### Product images:



[View online >](#)