

Product datasheet for **GA105884**

Human CD133 (PROM1) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PROM1 CRISPRa kit - CRISPR gene activation of human prominin 1 |
| Symbol: | PROM1 |
| UniProt ID: | <u>O43490</u> |
| Synonyms: | AC133; CD133; CORD12; MCDR2; MSTP061; PROML1; RP41; STGD4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)