

## Product datasheet for **GA105885**

### Human HM74 (HCAR3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HCAR3 CRISPRa kit - CRISPR gene activation of human hydroxycarboxylic acid receptor 3
Symbol:	HCAR3
UniProt ID:	<a href="#">P49019</a>
Synonyms:	GPR109B; HCA3; HM74; Puma-g; PUMAG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)