

Product datasheet for **GA105888**

Human TSC22D1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSC22D1 CRISPRa kit - CRISPR gene activation of human TSC22 domain family member 1
Symbol:	TSC22D1
UniProt ID:	Q15714
Synonyms:	Ptg-2; TGFB114; TSC22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)