

Product datasheet for **GA105890**

Human p35 (CDK5R1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CDK5R1 CRISPRa kit - CRISPR gene activation of human cyclin dependent kinase 5 regulatory subunit 1
Symbol:	CDK5R1
UniProt ID:	Q15078
Synonyms:	CDK5P35; CDK5R; NCK5A; p23; p25; p35; p35nck5a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)