

Product datasheet for **GA105892**

Human **PAG3 (ASAP2)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	ASAP2 CRISPRa kit - CRISPR gene activation of human ArfGAP with SH3 domain, ankyrin repeat and PH domain 2
Symbol:	ASAP2
UniProt ID:	O43150
Synonyms:	AMAP2; CENTB3; DDEF2; PAG3; PAP; Pap-alpha; SHAG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)