

## Product datasheet for **GA105925**

### Human EIF2S2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | EIF2S2 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 2 subunit beta                                                                                                                                                                                                  |
| Symbol:       | EIF2S2                                                                                                                                                                                                                                                                                                        |
| Components:   | GA105925G1, EIF2S2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA105925G2, EIF2S2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA105925G3, EIF2S2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>P20042</u>                                                                                                                                                                                                                                                                                                 |
| Synonyms:     | eIF-2-beta; EIF2; EIF2B; EIF2beta; PPP1R67                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)