

## Product datasheet for **GA105947**

### Human SNURF activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SNURF CRISPRa kit - CRISPR gene activation of human SNRPN upstream reading frame
Symbol:	SNURF
UniProt ID:	<u><a href="#">Q9Y675</a></u>
Synonyms:	SNRPN upstream reading frame; SNRPN upstream reading frame protein
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)