

Product datasheet for **GA105959**

Human TOP3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TOP3B CRISPRa kit - CRISPR gene activation of human DNA topoisomerase III beta
Symbol:	TOP3B
UniProt ID:	<u>O95985</u>
Synonyms:	TOP3B1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)