

Product datasheet for **GA105963**

Human Beta TRCP (BTRC) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | BTRC CRISPRa kit - CRISPR gene activation of human beta-transducin repeat containing E3 ubiquitin protein ligase |
| Symbol: | BTRC |
| UniProt ID: | Q9Y297 |
| Synonyms: | BETA-TRCP; betaTrCP; bTrCP; bTrCP1; FBW1A; FBXW1; FBXW1A; FWD1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)