

Product datasheet for **GA105964**

Human H3FI (HIST1H3F) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIST1H3F CRISPRa kit - CRISPR gene activation of human histone cluster 1 H3 family member f
Symbol:	H3C7
UniProt ID:	P68431
Synonyms:	H3/i; H3FI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)