

Product datasheet for **GA105976**

Human HSPB3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HSPB3 CRISPRa kit - CRISPR gene activation of human heat shock protein family B (small) member 3
Symbol:	HSPB3
Components:	GA105976G1, HSPB3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105976G2, HSPB3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105976G3, HSPB3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q12988
Synonyms:	DHMN2C; HMN2C; HSPL27
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)