

Product datasheet for **GA105987**

Human Proteinase Activated Receptor 4 (F2RL3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	F2RL3 CRISPRa kit - CRISPR gene activation of human F2R like thrombin or trypsin receptor 3
Symbol:	F2RL3
UniProt ID:	<u>Q96RI0</u>
Synonyms:	PAR4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)