

Product datasheet for **GA105988**

Human TAF1C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF1C CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor, RNA polymerase I subunit C
Symbol:	TAF1C
UniProt ID:	Q15572
Synonyms:	MGC:39976; SL1; TAFI95; TAFI110
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)