

## Product datasheet for **GA105989**

### Human TAF1B activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAF1B CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor, RNA polymerase I subunit B
Symbol:	TAF1B
UniProt ID:	<a href="#">Q53T94</a>
Synonyms:	MGC:9349; RAF1B; RAFI63; SL1; TAFI63
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)