

Product datasheet for **GA105990**

Human TAF1 48 (TAF1A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF1A CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor, RNA polymerase I subunit A
Symbol:	TAF1A
UniProt ID:	Q15573
Synonyms:	MGC:17061; RAFI48; SL1; TAFI48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)