

Product datasheet for **GA105999**

Human HIP1R activation kit by CRISPRa

Product data:

| | |
|---------------|----------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | HIP1R CRISPRa kit - CRISPR gene activation of human huntingtin interacting protein 1 related |
| Symbol: | HIP1R |
| UniProt ID: | <u>O75146</u> |
| Synonyms: | HIP3; HIP12; ILWEQ |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)