


Product datasheet for **GA106003**

Human **TM4SF5** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TM4SF5 CRISPRa kit - CRISPR gene activation of human transmembrane 4 L six family member 5
Symbol:	TM4SF5
UniProt ID:	<u>O14894</u>
Synonyms:	tetraspan transmembrane protein L6H; transmembrane 4 L six family member 5; transmembrane 4 superfamily member 5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)