

Product datasheet for **GA106004**

Human PKD2L1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PKD2L1 CRISPRa kit - CRISPR gene activation of human polycystin 2 like 1, transient receptor potential cation channel |
| Symbol: | PKD2L1 |
| UniProt ID: | Q9P0L9 |
| Synonyms: | PCL; PKD2L; PKDL; TRPP3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)