

Product datasheet for **GA106007**

Human PNR (TAAR5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAAR5 CRISPRa kit - CRISPR gene activation of human trace amine associated receptor 5
Symbol:	TAAR5
UniProt ID:	<u>O14804</u>
Synonyms:	PNR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)