


Product datasheet for **GA106011**

Human TAF172 (BTAF1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	BTAF1 CRISPRa kit - CRISPR gene activation of human B-TFIID TATA-box binding protein associated factor 1
Symbol:	BTAF1
UniProt ID:	O14981
Synonyms:	MOT1; TAF(II)170; TAF172; TAFII170
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)