

## Product datasheet for **GA106018**

### Human **PSTPIP1** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSTPIP1 CRISPRa kit - CRISPR gene activation of human proline-serine-threonine phosphatase interacting protein 1
Symbol:	PSTPIP1
UniProt ID:	<u><a href="#">O43586</a></u>
Synonyms:	CD2BP1; CD2BP1L; CD2BP1S; H-PIP; PAPAS; PSTPIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)