

## Product datasheet for **GA106019**

### Human RAI3 (GPCR5A) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	GPCR5A CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor class C group 5 member A
Symbol:	GPCR5A
UniProt ID:	<a href="#">Q8NFJ5</a>
Synonyms:	GPCR5A; PEIG-1; RAI3; RAIG1; TIG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)