

Product datasheet for **GA106043**

Human XKRY activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	XKRY CRISPRa kit - CRISPR gene activation of human XK related, Y-linked
Symbol:	XKRY
Components:	GA106043G1, XKRY gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106043G2, XKRY gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106043G3, XKRY gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2RUG3</u>
Synonyms:	XKRY1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

