

## Product datasheet for **GA106050**

### Human PIGQ activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIGQ CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class Q
Symbol:	PIGQ
UniProt ID:	<u><a href="#">Q9BRB3</a></u>
Synonyms:	c407A10.1; GPI1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

