

## Product datasheet for **GA106053**

### Human Retinal protein 4 (UNC119) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | UNC119 CRISPRa kit - CRISPR gene activation of human unc-119 lipid binding chaperone |
| Symbol:       | UNC119   |
| UniProt ID:   | <a href="#">Q13432</a>   |
| Synonyms:     | HRG4; IMD13; POC7; POC7A   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)              |

#### Product images:



[View online >](#)