

## Product datasheet for **GA106085**

### Human **KCNQ4** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | KCNQ4 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily Q member 4 |
| Symbol:       | KCNQ4  |
| UniProt ID:   | <a href="#">P56696</a>   |
| Synonyms:     | DFNA2; DFNA2A; KV7.4   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                  |

#### Product images:



[View online >](#)