

Product datasheet for **GA106094**

Human Synaptogyrin 3 (SYNGR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYNGR3 CRISPRa kit - CRISPR gene activation of human synaptogyrin 3
Symbol:	SYNGR3
UniProt ID:	O43761
Synonyms:	MGC20003
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)