

Product datasheet for **GA106107**

Human PCSK7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCSK7 CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 7
Symbol:	PCSK7
UniProt ID:	Q16549
Synonyms:	LPC; PC7; PC8; SPC7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)