

Product datasheet for **GA106115**

Human ST2 (IL1RL1) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | IL1RL1 CRISPRa kit - CRISPR gene activation of human interleukin 1 receptor like 1 |
| Symbol: | IL1RL1 |
| UniProt ID: | Q01638 |
| Synonyms: | DER4; FIT-1; IL33R; ST2; ST2L; ST2V; T1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)