

## Product datasheet for **GA106117**

### Human HTR3B activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | HTR3B CRISPRa kit - CRISPR gene activation of human 5-hydroxytryptamine receptor 3B |
| Symbol:       | HTR3B   |
| UniProt ID:   | <u><a href="#">O95264</a></u>   |
| Synonyms:     | 5-HT3B  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)             |

#### Product images:



[View online >](#)